



938 South Andreasen Drive, Suite E  
Escondido, CA 92029 USA

Phone: +1-760-738-1508  
Fax: +1-760-738-1530  
info@auroramedbiochem.com  
[www.auroramedbiochem.com](http://www.auroramedbiochem.com)

## EasySTOP™ Phosphatase Inhibitor Cocktail Set III

**ALTERNATE NAME:** Phosphatase Inhibitor Cocktail

**CATALOG #:** AUR10787

### Composition of Cocktail

Phosphatase Inhibitor	Mol. Wt.	Concentration In the Vial	Target Phosphatase
(-)-p-Bromotetramisole Oxalate	373.2	2.5 mM	Alkaline Phosphatase
Cantharidin	196.2	500 µM	Serine/Threonine Phosphatase
Calyculin A	1009.2	1 µM	Serine/Threonine Phosphatase

**FORM:** Liquid

**STORAGE:** Store at -20°C; Protect from air and moisture

**DESCRIPTION:** EasySTOP™ Phosphatase Inhibitor Cocktail Set III is a cocktail of three phosphatase inhibitors (see the table above) formulated in DMSO for the inhibition of both serine/threonine and alkaline phosphatases. Available as a 1 ml vial or as a set of five 1 ml vials. Each vial contains a 1 ml solution with the following components: 500 µM Cantharidin, 2.5 mM (-)-p-Bromotetramisole oxalate, and 1 µM Calyculin A, *Discodermia calyx*.

**HANDLING:** Do not take internally. Wear gloves and mask when handling the product! Avoid contact by all modes of exposure